

Product Data Sheet

Anti-Horse IgE (epsilon chain specific) Antibodies and conjugates

Cat# 80870	Goat Anti-Horse IgE, unlabeled	Size: 0.1 mg	Form: Liquid	Powder
Cat# 80871	Goat Anti-Horse IgE-HRP conjugate	Size: 0.5 ml	Form: Liquid	Powder

Description

Goats were immunized with antigen grade **Horse IgE**. Antibodies have been isolated using ammonium sulfate, ion-exchange, and affinity chromatography to remove non-specific antibodies. Specificity of antibody has been tested using IEP, immunodiffusion, and ELISA. The product reacts with Horse IgE with no significant reactivity with Horse IgG, or other serum proteins. The antibody may recognize other species IgE. Purified anti-Horse IgE antibody is supplied as unlabeled or conjugated to Biotin, FITC, AP or HRP.

Form and Storage

Cat# 10159, unlabeled

The antibody is supplied in PBS, pH 7.4, and 0.05% azide in either **lyophilized** (0.1 mg) or **liquid** form (0.1 mg/0.1 ml or or see lot sp. conc on the vial). Reconstitute powder in PBS to prepare stock solution. Store at -20oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Cat# 10160, HRP-conjugate

Purified antibody was coupled to HRP (RZ>3.0) using periodate method. The molar enzyme to protein (E/P) ratio = 4.0. The antibody is supplied in stabilizing buffer, 0.1% prolcin-300 as preservative in either **lyophilized** (0.5 ml) or **liquid** form. Reconstitute powder in PBS in 0.5 ml. Store at 4oC in suitable aliquots. Stability is ~6-12 months. Do not freeze and thaw.

Suggested conjugate dilutions are 1:1,000-1:10,000 ELISA, 1:1K-1:5K for western, and 1:200-1:1000 (IHC).

Recommended Working Dilution for ELISA

Working dilution for the specific application should be determined by the investigator to obtain the best conditions. Working solution should be prepared immediately before use and diluted solution should be discarded.

Related Material available for ADI

Horse IgE, IgG, IgE, and IgE ELISA kits

Anti-Horse IgE, IgG, IgM, IgA antibodies and conjugates

ELISA kits for the detection of Horse Antibodies

Western blot recycling kit (Use the same blot/strip to probe with multiple antibodies) **Strip Antibodies in 5 minutes** (no heating or mercaptoethanol)

MSDS: This product is considered non-hazardous as defined by OSHA (CFR 1910.1200. Nov 25, 1983). It can be disposed of in the drain.

80870-71

1200619A