

Product Specification Sheet

Rat Alpha 1 acid glycoprotein (A1AG) Protein

Cat. # A1AG19-N-25

Rat Alpha 1 acid glycoprotein (A1AG) Protein, pure

SIZE: 25 ug

A1AG serum protein, also called orosomucoid, is a monomer about 210 amino acid residues long; the amino acid sequence has been determined through 192 amino acids. A1AG is found in Cat plasma in concentrations of 55-140 mg per 100 ml. It is a classical standard glycoprotein (40 Kda, ~40% carbohydrates) for studies on the structure of the oligosaccharide units. Its biological significance is unknown, although it can bind progesterone 15 times as strongly as albumin. Sialic-acid-deficient alpha-1-AG has an affinity for vitamin B-12. Clinically, alpha 1 acid glycoprotein is an acute-phase reactant that together with haptoglobin is an indicator of acute inflammation. The alpha 1 acid glycoprotein:haptoglobin ratio is useful in studies of bone marrow disorders, hemolytic processes and metastases.

It is a positive acute phase reactant, the levels of which increase significantly in cats in response to inflammation, tissue injury and disease. A1AGP is purified using a combination of acid extraction, precipitation, gel filtration and affinity chromatography steps.

OROSOMUCOID 1; ORM1

Alternative titles; symbols

ORM

GLYCOPROTEIN, ALPHA-1-ACID, OF SERUM

ALPHA-1-ACID GLYCOPROTEIN

ALPHA-1-AGP; AGP1

Gene map locus 9q34.1-q34.3

Source of Antigen and Antibodies

A1AG protein (**cat# A1AG19-N-25**) was purified (>95%, ~50 Kda) from rat plasma. Purified A1AG protein is formulated in 150mM Tris, pH 7.4, and 0.05% azide (see lot specific protein concn on the vial). Lyophilized protein should be dissolved in PBS at a desired protein concentration. Store liquid or lyophilized protein at –20oC in suitable size aliquots. Do not freeze and thaw.

Storage

Short-term: unopened, undiluted vials for less than a week at 4oC.

Long-term: at –20C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

Stability: 6-12 months at –20oC or below.

Shipping: 4oC for solutions and room temp for lyophilized items.

General References: Dente L (1985) Nucl. Acid. Res. 13, 3941; Dente L (1987) EMBO J. 6, 2289; Board PG (1986) Gene 44, 127-131; Tokita K (1963) Nature, 200,266

*This product is for in vitro research use only.

Some New Antibodies from ADI...

Catalog#	ProdDescription
A1AG11-A	Anti-Cat Alpha 1-Acid Glycoprotein (A1-AGP) protein IgG, unlabeled
A1AG11-BT	Anti-Cat Alpha 1-Acid Glycoprotein (A1-AGP) protein IgG, Biotin conjugate
A1AG11-C	Cat Alpha 1-Acid Glycoprotein (A1-AGP) protein control for WB
A1AG11-HRP	Anti-Cat Alpha 1-Acid Glycoprotein (A1-AGP) protein IgG, HRP conjugate
A1AG19-N	Alpha 1 Acid Glycoprotein (A1-AGP), Cat Plasma
A1AG16-N-25	Alpha-1-Acid Glycoprotein Protein, Cat, purified
A1AG17-A	Anti-Dog Alpha-1-acid Glycoprotein IgG, aff pure
A1AG17-C	Alpha-1-Acid Glycoprotein Protein, Dog, control for western
A1AG17-N-25	Alpha-1-Acid Glycoprotein Protein, Dog, purified
A1AG18-A	Anti-Mouse a-1-Acid Glycoprotein IgG, aff pure
A1AG18-C	Alpha-1-Acid Glycoprotein Protein, Mouse, control for western
A1AG18-N-25	Alpha-1-Acid Glycoprotein Protein, Mouse, Purified
A1AG19-A	Anti-Rat a-1-Acid Glycoprotein IgG, aff pure
A1AG19-C	Alpha-1-Acid Glycoprotein Protein, Rat, purified control for western
A1AG19-N-25	Alpha-1-Acid Glycoprotein Protein, Rat, purified

ELISA kits

610-410-MLC Rat a-1-Acid Glycoprotein ELISA kit

A1AG19-N-25

110107A

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com