

Product Data Sheet

Cat # AB-10110

Mouse Anti HA-tag peptide IgG

Size: 100ug

Type:

Mouse Antibody Monoclonal

Immunogen:

HA-tag conjugated proteins

Ig Subclass:

Mouse IgG2b/kappa; purification by Protein A & ion exchange column; 1mg/mL in PBS after reconstitution; Reconstitute with H2O to get a 1 mg/ml concentration. Mix gently, wash the sides of the vial and wait 30-60 seconds before use.

Storage Procedures:

In lyophilized form, for long periods, store at 40C in a dry environment. After reconstitution, if not intended for use within a month, aliquot and store at -20C. Two years lyophilized, one month in solution at 40C. Antibody is shipped lyophilized at ambient temperature.

Titer:

By direct ELISA, 1 : 20,000 dilution will yield >1 O.D using alkaline phosphatase conjugated rabbit anti-mouse Ig (Jackson Laboratories).

Usage:

This item is for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals. If supplied in powder then reconstitute it in 100 ul water for 1 mg/ml stock and store in liquid at 40C for ~1 week or aliquots in suitable size and store at -20oC for long term storage..

Rev. 101001V

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com