

Product Data Sheet


---

Cat # ABT-400-100	Isooctanoate, Sodium (Pharma Grade)	Size: 100 g
Cat # ABT-400-500	Isooctanoate, Sodium (Pharma Grade)	Size: 500 g

---

General Information

Sodium Isooctanoate is known for its drying properties and good storage stability. It is used in the production of some coatings, as well as an accelerant for unsaturated polyester resins. It is also a perfect material for ampicillin and carbenicillin.



Molecular Formula:  
(C<sub>7</sub>H<sub>15</sub>COO)Na  
Molecular Weight:  
g/mol

**ABT-400-100**

**110520AO**

---

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)