

Sample Extraction Bag, SEB4



Part Number: ACC 00958/0025

Overview

Intended Use:

These sample mesh bags come filled with SEB4 sample extraction buffer. The bag provides an easy and effective way to grind plant samples. Simply cut open the bag, add your sample, and rub the outside of the bag with a pen or roller.

The special mesh lining of the bag aids in the maceration of plant tissues and the heavy plastic prevents the bag from ripping.

Includes:

- 25 Sample bags, mesh containing SEB4, 3 mL

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com