

**Agouti Related Protein (AGRP) Antibodies**

<b>Cat #</b> AGRP11-S	Rabbit Anti-Mouse AGRP, antiserum	<b>SIZE:</b> 100 ul
<b>Cat #</b> AGRP11-A	Rabbit Anti-Mouse AGRP IgG (aff pure)	<b>SIZE:</b> 100 ug
<b>Cat #</b> AGRP11-P	Mouse AGRP11 Control/blocking peptide	<b>SIZE:</b> 100 ug

Several common diseases such as type II diabetes, hypertension, cardiovascular diseases, hyperlipidemia, and some cancers are associated with obesity. An abnormal increase in body fat relative to lean tissue mass has been used as an indicator of obesity. High fat diet, certain environmental factors, and genetic linkage are the primary causes of obesity. In order to understand the genetic basis of obesity, several monogenic murine obesity models have been characterized including obese (*Ob*), diabetes (*db*), fat (*fat*), agouti yellow (*A<sup>y</sup>*), and tubby (*tub*). More recently, *Tub*, the human homolog of mouse *Tub*, *TULP1* & *TULP2* (for Tubby Like Proteins) and Agouti related protein (*AGRP*) have been cloned. The obesity associated with *Ay* mice may be due to ectopic expression of a secreted protein Agouti. *Agouti* protein (132 aa in human) is normally expressed in skin but its ubiquitous expression causes obesity. Agouti is a paracrine-signaling molecule that affects pigmentation by inhibiting the melanocortin receptor 1 (*MCR-1*). However, recombinant Agouti protein also antagonizes the *MC2R* and *MC4R*. *AGRP* (132 aa in human, chromosome 16q21) is normally expressed in adrenal and hypothalamus. *AGRP* levels are increased several folds in *ob/ob* mice. *AGRP* is a strong antagonist of *MC3R* and *MC4R*. Ubiquitous expression of *AGRP* in transgenic mice causes obesity without altering skin pigmentation.

**Source of Antigen and Antibodies**

<b>Antigen</b>	15-aa peptide of <b>Mouse AGRP (1); Designated (AGRP11-P or control peptide) /blocking peptide</b> conjugated to KLH; Epitope location ~C-terminus
<b>Ab Host/type</b>	Rabbit, Polyclonal unpurified antiserum ( <b>#AGRP11-S</b> ) and IgG, purified over antigen-agarose (Cat # <b>AGRP11-A</b> )
<b>2-Ab</b>	Cat # 20320, <b>goat anti-rabbit IgG-HRP</b> (AP, biotin, FITC conjugates also available).
<b>-ve control</b>	# 20009-1, Rabbit (non-immune) IgG, purified, suitable for ELISA, Western, IHC as -ve control

**Form & Storage of Antibodies/Peptide Control**

**Antiserum (unpurified)**

100ul solution lyophilized powder  
Supplied in Buffer: 0.05% azide  
**Reconstitute powder in 100 ul PBS**

**Affinity pure IgG**

100 ug/100ul solution lyophilized powder  
Supplied in **Buffer:** PBS+0.1% BSA  
**Reconstitute powder in PBS at 1mg/ml**

**Control/blocking peptide**

100 ug/100 ul solution lyophilized powder  
Supplied in Buffer: PBS pH 7.5,  
**Reconstitute powder in PBS at 1 mg/ml.**

**Storage**

**Short-term:** unopened, undiluted liquid vials at 20°C and powder at 4°C or -20°C..

**Long-term:** at -20°C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

**Stability:** 6-12 months at -20°C or below.

**Shipping:** 4°C for solutions and room temp for powder

**Recommended Usage**

**Western Blotting** (1:1K-5K for neat serum and 1-10 ug/ml for affinity pure using Chemiluminescence technique).

**ELISA** (1:10K-1:100K; using 50-100 ng of control peptide/well).

**Histochemistry & Immunofluorescence:** not tested. We recommend the use of affinity pure antibody at 2-20 ug/ml.

**Specificity & Cross-reactivity**

Mouse AGRP11-P peptide has homology with human and bovine (85%). No homology is seen with Agouti, tubby, *tub*, *TULP1*, and *TULP2*. Antibody crossreactivity in various species is not known. The **AGRP11-P**, Control peptide, because of its low mol. Wt (<3 kDa), is not suitable for Western. It should be used for ELISA or antibody blocking experiments (use 5-10 ug control peptide per 1 ug of aff pure IgG or 1 ul antiserum) to confirm antibody specificity.

**General References:** 1. Shutter RJ et al (1997) *Genes Develop.* 11, 593; Ollmann MM et al (1997) *Science* 278, 135.

(2) Citations of ADI's Antibodies (see web site for updated list)

Bertorelli R, 2005, *Behavioural Brain Research*, 167, 55-62, IHC

*\*This product is for In vitro research use only.*

**Related material available from ADI**

Anti-Agouti, AGRP, Tubby, *TUB*, *TULP1*, *TULP2*, Leptin, and Melanocortin receptors 91-5)

AGRP11-S-A-P 71129S

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)