

Product Specification Sheet

**Amphetamine Antibody**

|  |   |
|--|---|
| <b>Cat. #</b>                          | AMPT11-M  |
| <b>Description</b>                     | Mouse Monoclonal Anti-Amphetamine IgG (aff pure)  |
| <b>Host Animal</b>                     | Mouse   |
| <b>Isotype</b>                         | IgG2b   |
| <b>Form</b>                            | Protein A/G affinity purified from ascites fluid supplied in PBS, pH 7.4 and 0.05% sodium azide as preservative   |
| <b>Antibody concentration</b>          | 1 mg/ml   |
| <b>Storage</b>                         | Short term at 2-8oC; Long term at -20oC   |
| <b>Stability</b>                       | ~1 year or more at -20oC; Avoid repeated freeze and thaw.   |
| <b>Specificity and crossreactivity</b> | Specifically Reacts with amphetamine  |
| <b>Suggested Usage</b>                 | Suitable for ELISA as coating antibody; concentrations and coating conditions must be optimized for the assay.    |
| <b>MSDS Info</b>                       | Contains minimal concn of sodium azide (0.05%). Consult the handling and disposal recommendations of sodium azide |

AMPT11-M

70314A

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)