

Product Specification Sheet

Amphetamine Antibody

Cat. #	AMPT11-M
Description	Mouse Monoclonal Anti-Amphetamine IgG (aff pure)
Host Animal	Mouse
Isotype	IgG2b
Form	Protein A/G affinity purified from ascites fluid supplied in PBS, pH 7.4 and 0.05% sodium azide as preservative
Antibody concentration	1 mg/ml
Storage	Short term at 2-8oC; Long term at -20oC
Stability	~1 year or more at -20oC; Avoid repeated freeze and thaw.
Specificity and crossreactivity	Specifically Reacts with amphetamine
Suggested Usage	Suitable for ELISA as coating antibody; concentrations and coating conditions must be optimized for the assay.
MSDS Info	Contains minimal concn of sodium azide (0.05%). Consult the handling and disposal recommendations of sodium azide

AMPT11-M

70314A

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com