

Product Specification Sheet

**Angiopoietin-1 (Ang-1) Antibodies**

<b>Cat. ANG11-S</b>	Rabbit Anti-Mouse Ang-1 Antiserum # 1	<b>SIZE:</b> 100 ul
<b>Cat. ANG11-A</b>	Rabbit Anti-Mouse Ang-1 IgG # 1 (aff pure)	<b>SIZE:</b> 100 ug
<b>Cat. ANG11-P</b>	Mouse Ang-1 Control/blocking peptide # 1	<b>SIZE:</b> 100 ug

Embryonic vascular system undergoes a series of complex, highly regulated series of events involving differentiation, migration and association of primitive endothelial cells. This process is termed vasculogenesis. A further remodeling of the primitive vascular system forms the mature cardiovascular system. This process is known as **angiogenesis** (sprouting of new capillary vessels from pre-existing vasculature). Angiogenesis accounts for the formation of vasculature into previously avascular organs such as brain and kidney. Angiogenic activity in the adult is required during the normal tissue repair, and for the remodeling of the female reproductive organs (ovulation and placental development).

**Angiopoietin-1** (mouse and human **Ang-1**; 498 AA; ~ 98% identity) is an angiogenic secreted protein that interact with endothelial specific Tie-2 receptor. It is primarily expressed in developing endothelial cells. During embryonic development, Ang-1 binds and induces tyrosine phosphorylation of Tie-2. Ang-1 appears to play a crucial role in mediating matrix and mesenchyme. It mediates blood vessel maturation and stability. Ang-1 may play a critical role in heart development. Ang-1 deficient mice mimic the phenotype exhibited by animals deficient in Tie-2. A homolog of Ang-1, termed **Angiopoietin-2** (mouse and human **Ang-2**, 496 AA; ~85% identity) has recently been identified. It may act an antagonist for Ang-1 and Tie-2. Ang-1 and Ang-2 have ~60% sequence homology.

**Source of Antigen and Antibodies**

<b>Antigen</b>	20-aa peptide from mouse Ang-1 (1) KLH conjugate ( <b>Cat. ANG11-P</b> ). Epitope location ~ N-terminus
<b>Ab Host/type</b>	Rabbit, polyclonal Unpurified antiserum <b>Cat # ANG11-S</b> and Aff pure IgG ( <b>cat #ANG11-A</b> )
<b>2ab</b>	Cat # 20320, goat anti-rabbit IgG-HRP (AP, biotin, FITC conjugates also available
<b>-ve control</b>	# 20009-1, Rabbit (non-immune) IgG, purified, suitable for ELISA, Western, IHC as -ve control

**Form & Storage of Antibodies/Peptide Control**

**Antiserum (unpurified)**  
100ul solution lyophilized powder  
Supplied in Buffer: 0.05% azide  
**Reconstitute** powder in 100 ul PBS

**Affinity pure IgG**  
100 ug/100ul solution lyophilized powder  
Supplied in **Buffer:** PBS+0.1% BSA  
**Reconstitute** powder in PBS at 1mg/ml

**Control/blocking peptide**  
100 ug/100 ul solution lyophilized powder  
Supplied in Buffer: PBS pH 7.5,  
**Reconstitute** powder in PBS at 1 mg/ml.

**Storage**

**Short-term:** unopened, undiluted liquid vials at 20°C and powder at 4°C or -20°C..

**Long-term:** at -20°C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

**Stability:** 6-12 months at -20°C or below.

**Shipping:** 4°C for solutions and room temp for powder

**Recommended Usage**

**Western Blotting** (1:1K-5K for neat serum and 1-10 ug/ml for affinity pure using Chemiluminescence technique). Expected size of Ang-1 is ~55 kDa. An antibody made to Ang-1 epitope has detected ~ 70 kDa glycosylated, secreted protein in the conditioned media from cells expressing recombinant Ang-1 (1). (refs 2)

**ELISA** (1:10K-1:100K; using 50-100 ng of control peptide/well).

**Histochemistry & Immunofluorescence:** We recommend the use of affinity pure antibody at 2-20 ug/ml (refs 2).

**Specificity & Cross-reactivity**

The ANG11-P control immunogenic peptide sequences is 90% conserved in human, 95% in pig, rat and bovine Ang-1. Ang-1 control peptide has no significant sequence homology with Ang-2. Antibody cross-reactivity in various species has not been studied. Control peptide, because of its low mol. Wt (<3 kDa), is not suitable for Western. It should be used for ELISA or antibody blocking experiments (use 5-10 ug control peptide per 1 ug of aff pure IgG or 1 ul antiserum) to confirm antibody specificity (see detailed protocol at: the web site).

**General References:** (1) Davis S (1996) Cell 87, 1161; Nomura N (1994) DNA Res. 1, 27-35; Suri, S (1996) Cell 87, 1171; Davis, S (1996) cell 87, 1161; Maisonpierre C et al (1997) Science 277, 55.

**(2) Citations of ADI's Antibodies** (see web site for updated list)

Yuan Hai Tao, 2002 Kidney Intl. 61:2078, WB, IHC  
Long, D A, 2001, J Am Soc Nephrol 12: 2721, WB, IHC  
Zheng W, 2004, Am J Physiol Heart Circ Physiol. 286: 1994, WB  
Krikun G, 2002, Am. j. Pathol. 161, 979, WB, IHC  
Hayashi, T. 2003, J. Cerebral Blood Flow & Metabol. 23:166, IHC  
Long, DA, 2001, J Am Soc Nephrol 2001 12: 2721, WB, IHC

*\*This product is for in vitro research use only.*

**Related material available from ADI**

ANG11-S-A-P

71216S

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)