



Viale Piceno 137/f
61032 Fano PU (IT)
Telephone + 39 (0) 721830605
FAX +39 (0) 721837154
e-mail: info@diatheva.com
www.diatheva.com

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura,
Delhi – 110034 (INDIA).
Ph: +91-11-42208000, 42208111, 42208222
Mobile: +91-9810521400
Fax: +91-11-42208444
Email: customerservice@lifetechindia.com
Web: www.lifetechindia.com

Rabbit Anti Listeriolysin O (LLO)

ANT0006
200µl

Description	Listeriolysin O (LLO), a major virulence factor produced by all pathogenic strains of <i>L. monocytogenes</i> , has been identified as a candidate antigen for a serological assay.
Product type	Polyclonal antibody
Immunogen	Recombinant Listeriolysin O protein expressed in <i>E. coli</i>
Source	Rabbit
Reacts with	Wild type Listeriolysin O protein from <i>Listeria monocytogenes</i>
Specificity	NA
Tested applications	Western Blotting, ELISA
Recommended dilutions	<p>Recommended starting dilutions can vary lot-to-lot. Consult the product information label in the package for lot specific values.</p> <p>Note: When using any primary antibody or fluorescence-labelled secondary antibody for the first time, titrate out the antibody to determine which dilution allows the strongest specific signal with the lowest background for your sample.</p> <p><i>For untested applications or species please refer to the S.M.A.K. program.</i></p>
Purity	Polyclonal immunoglobulins purified by protein A affinity chromatography.
Form	Liquid. Supplied in 100mM sodium citrate, 50mM Tris and 0.05% v/v glycerol. Neutral pH.
Reconstitution and Storage	<p>Shipped at +4°C. When stored at +4°C, the antibody is stable for 18 months. For extended storage, the solution may be frozen at –20°C in working aliquots.</p> <p>Note: Avoid repeated freezing and thawing cycles.</p>
References	Available on library section: http://www.diatheva.com/library.htm