
Product Data Sheet

Cat#: AP-309-U
Cat#: AP-309-B, Biotin labeled
Cat#: AP-309-F, FITC labeled

Product Description: EphA2 Receptor (YSA peptide)

Aptamer Type: Peptide

Sequence: YSAYPDSVPMMS; 12-mer

Size 5 mg

Mol. Wt: 1483.71 g/mole

Purity: >95%

Affinity: 186 nM (reported value)

Comments: n/a

Notes: Peptide has ephrin-like activity that it stimulates EphA2 tyrosine phosphorylation and signaling. Furthermore, this peptide can deliver phage particles to endothelial and tumor cells expressing EphA2

References: Koolpe et al. "An Ephrin Mimetic Peptide That Selectively Targets the EphA2 Receptor." The Journal of Biological Chemistry, 277(2002): 46974-46979.

150813V