

Product Data Sheet

Cat#: AP-310-U
Cat#: AP-310-B, Biotin labeled
Cat#: AP-310-F, FITC labeled

Product Description: EphA2 Receptor (SWL peptide)

Aptamer Type: Peptide

Sequence: SWLAYPGAVSYR;12-mer

Size 5 mg

Mol. Wt: 1505.73 g/mole

Purity: >95%

Affinity: 678 nM (reported value)

Comments: n/a

Notes: Peptide has ephrin-like activity in that it stimulates EphA2 tyrosine phosphorylation and signaling. Furthermore, this peptide can deliver phage particles to endothelial and tumor cells expressing EphA2

References: Koolpe et al. "An Ephrin Mimetic Peptide That Selectively Targets the EphA2 Receptor." The Journal of Biological Chemistry, 277(2002): 46974-46979.

150813V

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com