

Product Data Sheet

Cat#: AP-321-U
Cat#: AP-321-B, Biotin labeled
Cat#: AP-321-F, FITC labeled

Product Description: Mouse Lung

Aptamer Type: Peptide

Sequence: CGFECVRQCPERC; 13-mer

Size 5 mg

Mol. Wt: 1684.03 g/mole

Purity: >95%

Affinity: n/a

Comments: n/a

Notes: The enrichment in homing to the target organs relative to an unselected phage was 3-35-fold. Peptide sequences that conferred selective phage homing to the vasculature of lung, skin, and pancreas were characterized in detail. Immunohistochemistry showed that the phage localized in the blood vessels of their target organ. When tested, the phage homing was blocked in the presence of the cognate peptide.

References: Rajotte et al. "Molecular Heterogeneity of the Vascular Endothelium Revealed by In Vivo Phage Display." *Journal of Clinical Investigations*, 102(1998): 430-437

150813V

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com