

Product Data Sheet

Cat#: AP-328-U
Cat#: AP-328-B, Biotin labeled
Cat#: AP-328-F, FITC labeled

Product Description: PC-3 (G1)

Aptamer Type: Peptide

Sequence: IAGLATPGWSHWLAL; 15-mer

Size 5 mg

Mol. Wt: 837.98 g/mole

Purity: >95%

Affinity: n/a

Comments: n/a

Notes: One resulting phage clone (G1) displaying the peptide sequence IAGLATPGWSHWLAL was fluorescently labeled with the near-infrared fluorophore AlexaFluor 680 and was evaluated both in vitro and in vivo for its ability to bind and target PC-3 prostate carcinomas. The fluorescently labeled phage clone (G1) had a tumor-to-muscle ratio of approximately 30 in experiments. In addition, prostate tumors (PC-3) were readily detectable by optical-imaging methods.

References: Beck et al. "Identification of a peptide that interacts with Nestin protein expressed in brain cancer cells." *Biomaterials*, 32(2011): 8518-8528.

150813V

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com