

Product Specification Sheet

Human Apolipoprotein C-II (Apo C-II) protein and antibodies

<input type="checkbox"/> Cat. # APOC22-A	Rabbit Anti-Human Apo C-II protein IgG, aff pure	SIZE: 100 ul
<input type="checkbox"/> Cat. # APOC22-C	Purified Human Apo C-II protein control for WB	SIZE: 100 ul

Apolipoproteins are proteins that bind to fats (lipids). They form lipoproteins, which transport dietary fats through the bloodstream. Dietary fats are digested in the intestine and carried to the liver. Fats are also synthesized in the liver itself. Fats are stored in fat cells (adipocytes). Fats are metabolized as needed for energy in the skeletal muscle, heart, and other organs and are secreted in breast milk. Apolipoproteins also serve as enzyme co-factors, receptor ligands, and lipid transfer carriers that regulate the metabolism of lipoproteins and their uptake in tissues.

Apolipoprotein C-II is a protein that in humans is encoded by the APOC2 gene. The protein encoded by this gene is secreted in plasma where it is a component of very low density lipoprotein. This protein activates the enzyme lipoprotein lipase, which hydrolyzes triglycerides and thus provides free fatty acids for cells. Mutations in this gene cause hyperlipoproteinemia type IB, characterized by hypertriglyceridemia, xanthomas, and increased risk of pancreatitis and early atherosclerosis. The concentration of Apo CII in normal plasma is 3-8 mg per 100 ml.

Source of Antigen, Antibodies, and positive controls

Antigen	Purified human Apo C-II protein
Ab Host/type	Rabbit, Poly IgG # APOC22-A
2-Ab	Goat Anti-rabbit IgG-HRP cat # 20320 (AP, biotin, FITC conjugates also available)
-ve control IgG	# 20009-1, Rabbit (non-immune) IgG, purified, suitable for ELISA, Western, IHC as -ve control

Cat# APOC22-C

For Western blot +ve control, purified human Apo C-II protein (Cat # **APOC22-C**) is supplied in SDS-PAGE sample buffer (reduced). Load 10 ul/lane of **APOC22-C** for good visibility with antibody Cat # **APOC22-A**. Store at -20oC in suitable size aliquots. SDS may crystallize in cold conditions. It should redissolve by warming before taking it from the stock. It should be heated once prior to loading on gels. If the product has been stored for several weeks, then it may be preferable to add 5 ul of fresh 2x sample buffer per 10 ul of the **APOC22-C** solution prior to

heating and loading on gels. This preparation is not biologically active. It is not suitable for ELISA or other applications where native protein is required. Do not freeze, thaw, or heat repeatedly.

Form & Storage of Antibodies/Peptide Control

IgG purified

100 ul/vial solution lyophilized powder

contains 0.05% sodium azide

Reconstitute powder 100 ul PBS

Recommended Usage

Western Blotting 1-3 ug/ml using Chemiluminescence technique). Human ApoC-II s approx. ~8.8 KDa (1).

ELISA (1:10K-1:100K; using 50-100 ng of control peptide/well).

Histochemistry & Immunofluorescence: not tested. We recommend a antibody testing at 5-20 ug/ml.

Specificity & Cross-reactivity

The antibody reacts with human Apo C-II in ELISA and Western. Antibody crossreactivity in various species is not established. **APOC22-C** protein control should be used a positive control for western.

General References: Jackson RI (1977) PNAS 74, 1942-1945; Hegley RA (1992) Dis. Markers 9, 73-80; Wei CF (1986) JBC 260, 15211-1521; Baggio G (1986) J. Clin. Invest. 77, 520-527.

*This product is for In vitro research use only.

Related material available from ADI

Ant-ApoA, ApoB, ApoC, ApoE -Beta amyloid 1-40, 1-42, APP, Parkin, Synucleins (α , β , γ), Presenilins 1, 2, ERAB

APOC22-A-C 100317A

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com