

Product Specification Sheet

**Aquaporin 5 (AQP5) Antibodies**

<b>Cat # AQP51-P</b>	Rat AQP5 Control/blocking peptide	<b>SIZE:</b> 100 ug
<b>Cat # AQP51-S</b>	Rabbit Anti-rat AQP5 antiserum #1	<b>SIZE:</b> 100 ul
<b>Cat # AQP51-A</b>	Rabbit Anti-rat AQP5 IgG #1 (aff pure)	<b>SIZE:</b> 100 ug

Water is a critical component of all living cells. Interestingly, tissue membranes show a great degree of water permeability. Mammalian red cells, renal proximal tubules, and descending thin limb of Henle are extraordinarily permeable to water. Water crosses hydrophobic plasma membranes either by simple diffusion or through a facilitative transport mechanism mediated by special protein "aquaporin". Over the last decade, genes for several members of aquaporin family (AQP1-9 and AQPAP) have been cloned, expressed, and their distribution studied in many tissues. **AQP5** (human, mouse, and rat 265 aa) is found in several epithelial tissues (salivary and lachrymal glands, corneal epithelium in the eye and lung). It is implicated in generation of saliva, tears, and pulmonary secretions. AQP families of proteins are predicted to contain six transmembrane domains.

**FUNCTION:** Forms a water-specific channel. Implicated in the generation of saliva, tears, and pulmonary secretions.

**SUBCELLULAR LOCATION:** Multi-pass membrane protein.

**TISSUE SPECIFICITY:** Salivary glands, lacrimal glands, corneal epithelium in eye, trachea and lung.

**DOMAIN:** Aquaporins contain two tandem repeats each containing three membrane-spanning domains and a pore-forming loop with the signature motif Asn-Pro-Ala (NPA).

**SIMILARITY:** Belongs to the MIP/aquaporin (TC 1.A.8)

Protein name Aquaporin-5 ; Synonym AQP-5

Gene name Name: Aqp5

**Source of Antigen and Antibodies**

<b>Antigen</b>	17-aa peptide of Rat AQP5 (Gene Accession #P47864; <b>Designated (AQP51-P or control peptide) conjugated to KLH; epitope location ~ C-terminus</b> )
<b>Ab Host/type</b>	Rabbit, polyclonal, Unpurified antiserum (cat # AQP51-S) Aff pure IgG1 ( <b>cat #AQP51-A</b> ) purified over antigen-agarose column
<b>2-Ab</b>	<b>Goat Anti-rabbit IgG-HRP</b> cat # 20320 (AP, biotin, FITC conjugates also available)
<b>-ve control</b>	<b># 20009-1, Rabbit (non-immune) IgG, purified, suitable for ELISA, Western, IHC as -ve control</b>

**Form & Storage of Antibodies/Peptide Control**

**Antiserum (unpurified)**

100ul solution lyophilized powder  
Supplied in Buffer: 0.05% azide  
**Reconstitute powder in 100 ul PBS**

**Affinity pure IgG**

100 ug/100ul solution lyophilized powder  
Supplied in **Buffer:** PBS+0.1% BSA  
**Reconstitute powder in PBS at 1mg/ml**

**Control/blocking peptide**

100 ug/100 ul solution lyophilized powder  
Supplied in Buffer: PBS pH 7.5,  
**Reconstitute powder in PBS at 1 mg/ml.**

**Storage**

**Short-term:** unopened, undiluted liquid vials at -20OC and powder at 4oC or -20oC..

**Long-term:** at -20C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

**Stability:** 6-12 months at -20oC or below.

**Recommended Usage**

We recommend the use of 0.5-1% milk in all primary/secondary antibody-enzyme conjugate incubations in order to suppress non-specific bands.

**Western blotting:** 1:1K-5K for neat serum and 1-10 ug/ml for affinity pure antibody using ECL (see published refs in 2).

**ELISA:** Control peptide can be used to coat ELISA plates at 1 ug/ml and detected with antibodies (1:10-50K for neat serum and 0.5-1 ug/ml for affinity pure).

**Histochemistry & Immunofluorescence:** We recommend the use of affinity-purified antibody at 2-10 ug/ml in paraformaldehyde fixed sections of tissues ((see published refs in 2).

**Specificity & Cross-reactivity**

The rat AQP51-P peptide is 100% conserved in mouse, and 80% in human, 76% in ovine AQP5. However, no significant homology is seen with other AQPs or any other protein. Antibody cross-reactivity in various species has not been studied. Control peptide, because of its low mol. Wt (<3 kDa), is not suitable for Western. It should be used for ELISA or antibody blocking experiments (use 5-10 ug control peptide per 1 ug of aff pure IgG or 1 ul antiserum) to confirm antibody specificity.

**General References:**

(1) Raina, S. et al. (1995) *J. Biol. Chem.* **270**:1908-1912; Lee MD et al (1996) *JBC* **271**, 8599-8604; Krane CM et al (1999) *Mamm. Genome* **10**, 498-505.

**Citations of ADI's antibodies for AQP5** (updated list at the web site).

Borok Z, 2002, *AJP Lung Cell Mol Physiol* **282**: L599, WB IF  
Blank M, 2004, *Mol. Cell. Neurosci.* **26**, 530-543, IF  
Chen J 2004, *Lab. Invest.* **84**, 727 - 735 WB  
Li J 2004, *Lab. Invest.* **84**, 1430 - 1438 WB IHC  
Borok Zea, 1998, *Am. J. Respir. Cell Mol. Bio.* **18**: 554-561 WB,  
Beroukas D, 2001, *Lancet* **358**, 1875-1876 IHC,  
Larina O 2005, *J. Cell Sci.* **118**: 4131 - 4139 IF  
Lindsay LA, 2004 *Acta Histochemica* **106**, 299-307, IF  
DaSilva N2005 *Biol Reprod.* Oct 2005 WB IF  
Rabb H, 2003 *Kidney Intl.* **63**(2):600-606 WB  
Wang W 2003 *Calcified T Intl.* **72**(3):222-7 WB

\*This product is for *In vitro* research use only.

AQP51-S-A-P 70910A

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)