

Product Specification Sheet

Classical Adjuvant Combo Pak-1

□ **Cat. #** AV-3000-PK-1 Classical Adjuvant Combo Pak-1 (contains 10 ml each CFA (#AV-3010-10); IFA (AV-3015-10), and Squalene (#AV-3020-10) **SIZE:** 1 PK

General Information: The word '**adjuvant**' is derived from the Latin word '**adjuvare**' which means '**to help**'. Therefore, Immunologic Adjuvants are added to vaccines to stimulate the immune system's response to the target antigen, but do not in themselves confer immunity. Adjuvants act in various ways in presenting an antigen to the immune system. Adjuvants can act as a depot for the antigen, presenting the antigen over a long period of time, thus maximizing the immune response before the body clears the antigen. Examples of depot type adjuvants are oil emulsions. Adjuvants can also act as an irritant which causes the body to recruit and amplify its immune response. A tetanus, diphtheria, and pertussis vaccine, for example, contains minute quantities of toxins/toxoids produced by each of the target bacteria. The body's immune system develops an antitoxin to the bacteria's toxins, not to the aluminum, but would not respond enough without the help of the aluminum adjuvant. Adjuvants have also evolved as substances that can aid in stabilizing formulations of antigens, especially for vaccines administered for animal health.

Adjuvants augment the effects of a vaccine by stimulating the immune system to respond to the vaccine more vigorously, and thus providing increased immunity to a particular disease. Adjuvants accomplish this task by mimicking specific sets of evolutionarily conserved molecules, so called PAMPs, which include liposomes, lipopolysaccharide (**LPS**), molecular cages for antigen, components of bacterial cell walls (e.g., **flagellins**), and endocytosed nucleic acids such as double-stranded RNA (**dsRNA**), single-stranded DNA (**ssDNA**), and unmethylated CpG dinucleotide-containing DNA (**ODNs**). Natural proteins such as **ovalbumin** or OVA-peptides and key hole limpet hemocyanins (**KLH**) are also being explored not only serve as carrier protein but also as adjuvants. Because immune systems have evolved to recognize these specific antigenic moieties, the presence of an adjuvant in conjunction with the vaccine can greatly increase the innate immune response to the antigen by augmenting the activities of dendritic cells (DCs), lymphocytes, and macrophages by mimicking a natural infection. Furthermore, because adjuvants are attenuated beyond any function of virulence, they pose little or no independent threat to a host organism.

For human vaccines, aluminum hydroxide (Alum) based adjuvants (Aluminum hydroxide or Alhydrogel; Aluminium phosphate or Adjuvax) are the only **FDA-approved adjuvants**. Vaccine components that are formulated in Alum are called "Adsorbed Vaccines". The effectiveness of each salt as an adjuvant depends on the characteristics of the specific vaccine and how the manufacturer prepares the vaccine

Not all vaccines contain Alum because an adjuvant may not have been needed, was not expected to increase the desired immune response, or was going to cause an imbalance in the immune response. For example, **inactivated Polio Virus (IPV/IPOL)** vaccine, measles, mumps and rubella vaccine (**MMR/MMRI/MMRV**), **Varicella or chickenpox vaccine (Varivax/Proquad/MMRV)**, **Meningococcal conjugate (MCV4/Menomune/Menactra)** vaccine, and **influenza vaccines (Fluzone/Flulaval/Flumist/Fluvirin etc)** do not contain aluminum salts.

Product Information

Freund's Adjuvants may be used to produce water-in-oil emulsions of immunogens. The mycobacteria in Complete Freund's adjuvant attract macrophages and other cells to the injection site which enhances the immune response. For this reason, Complete Freund's Adjuvant is used for initial injections and Incomplete Freund's Adjuvant for subsequent boosts

Squalene is a natural 30-carbon organic compound originally obtained for commercial purposes primarily from shark liver oil humans. Squalene as an adjuvant is manufactured by Novartis

Classical adjuvant combo pack Contains 3 most common classical oil-based adjuvant emulsions; 10 ml each of Complete Freund's Adjuvant/**CFA** (#AV-3010-10); Incomplete Freund's adjuvant/**IFA** (AV-3015-10), and **Squalene** (#AV-3020-10) to test adjuvantcy with a given antigen.

Designed to provide various formulations of classical oil based adjuvants for testing adjuvantcy with a given antigen.

Storage and Stability: Shipped at room temperature and it should be stored at room temp. DO NOT FREEZE. Stable for 6 months.

Suggested Usage: Adsorption of most proteins can be achieved by mixing equal volume of antigen and adjuvant. Protein left in the supernatant can be measured by protein assay or specific ELISA.

Related items :

Catalog#	ProdDescription
AV-3000-PK-1	Classical Adjuvant Combo Pak-1 (contains 10 ml each CFA (#AV-3010-10); IFA (AV-3015-10), and Squalene (#AV-3020-10)
AV-3010-10	Complete Freund's Adjuvant (CFA);vaccine adjuvant
AV-3010-100	Complete Freund's Adjuvant (CFA);vaccine adjuvant
AV-3015-10	Incomplete Freund's Adjuvant (IFA);vaccine adjuvant
AV-3015-100	Incomplete Freund's Adjuvant (IFA);vaccine adjuvant
AV-3020-10	Squalene (oil-in-water nano emulsion);vaccine adjuvant

Complete list is available at:
http://4adi.com/objects/catalog/product/extras/Vaccine_Adjuvants_flr.pdf

AV-3000-PK-1 140519P

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com