

Product Specification Sheet

Ovalbumin peptide OVA (323-339) MHC class II peptide

- Cat. #. AV-9320-1
 Cat. #. AV-9320-10

Ovalbumin peptide OVA (323-339) MHC class II peptide
Ovalbumin peptide OVA (323-339) MHC class II peptide

SIZE: 1 mg
SIZE: 10 mg

General Information: The word 'adjuvant' is derived from the Latin word 'adjuvare' which means 'to help'. Therefore, Immunologic Adjuvants are added to vaccines to stimulate the immune system's response to the target antigen, but do not in themselves confer immunity. Adjuvants act in various ways in presenting an antigen to the immune system. Adjuvants can act as a depot for the antigen, presenting the antigen over a long period of time, thus maximizing the immune response before the body clears the antigen. Examples of depot type adjuvants are oil emulsions. Adjuvants can also act as an irritant which causes the body to recruit and amplify its immune response. A tetanus, diphtheria, and pertussis vaccine, for example, contains minute quantities of toxins/toxoids produced by each of the target bacteria. The body's immune system develops an antitoxin to the bacteria's toxins, not to the aluminum, but would not respond enough without the help of the aluminum adjuvant. Adjuvants have also evolved as substances that can aid in stabilizing formulations of antigens, especially for vaccines administered for animal health.

Adjuvants augment the effects of a vaccine by stimulating the immune system to respond to the vaccine more vigorously, and thus providing increased immunity to a particular disease. Adjuvants accomplish this task by mimicking specific sets of evolutionarily conserved molecules, so called PAMPs, which include liposomes, lipopolysaccharide (LPS), molecular cages for antigen, components of bacterial cell walls (e.g., flagellins), and endocytosed nucleic acids such as double-stranded RNA (dsRNA), single-stranded DNA (ssDNA), and unmethylated CpG dinucleotide-containing DNA (ODNs). Natural proteins such as ovalbumin or OVA-peptides and key hole limpet hemocyanins (KLH) are also being explored not only serve as carrier protein but also as adjuvants. Because immune systems have evolved to recognize these specific antigenic moieties, the presence of an adjuvant in conjunction with the vaccine can greatly increase the innate immune response to the antigen by augmenting the activities of dendritic cells (DCs), lymphocytes, and macrophages by mimicking a natural infection. Furthermore, because adjuvants are attenuated beyond any function of virulence, they pose little or no independent threat to a host organism.

For human vaccines, aluminum hydroxide (Alum) based adjuvants (Aluminum hydroxide or Alhydrogel; Aluminium phosphate or Adjuvphos) are the only **FDA-approved adjuvants**. Vaccine components that are formulated in Alum are called "Adsorbed Vaccines". The effectiveness of each salt as an adjuvant depends on the characteristics of the specific vaccine and how the manufacturer prepares the vaccine

Not all vaccines contain Alum because an adjuvant may not have been needed, was not expected to increase the desired immune response, or was going to cause an imbalance in the immune response. For example, **inactivated Polio Virus** (IPV/IPOL) vaccine, measles, mumps and rubella vaccine (MMR/MMRI/MMRV), **Varicella or chickenpox vaccine** (Varivax/Proquad/MMRV), **Meningococcal conjugate** (MCV4/Menomune/Menactra) vaccine, and **influenza vaccines** (Fluzone/Flulaval/Flumist/Fluvirin etc) do not contain aluminum salts.

Product Information

Ovalbumin (abbreviated OVA) is the main protein found in egg white, making up 60-65% of the total protein. Ovalbumin displays sequence and three-dimensional homology to the serpin super family. The function of ovalbumin is unknown, although it is presumed to be a storage protein.

OVA 323-339 encompasses an allergenic and antigenic epitope of the ovalbumin protein1 .OVA 323-339 peptide binds to I-A(d) Major Histocompatibility Complex (MHC) class II protein2 . OVA 323-339 has been used extensively to study the nature of class II MHC-peptide binding and T-cell activation

Formula : C74H120N26O25

Mol.wt : 1773.9

Appearance: White to off-white powder

Solubility: 2 mg/ml in water

Storage and Stability: Shipped at room temperature and it should be stored at 2-8 C. Long term storage at -20 C for up to 6 months. Avoid repeated freeze thaw cycles.

References: Rotzschke O. et al., (1991). Eur J Immunol.21(11):2891-4. 2. Huntington JA; Stein PE (2001). *Journal of Chromatography B* 756 (1-2): 189-198

Related items:

Catalog#	ProdDescription
AV-9305-10	Keyhole Limpet Hemocyanin (Megathura Crenulata)
AV-9310-1000	Chicken egg ovalbumin protein (ELISA, antigen, allergy grade)
AV-9315-1	Ovalbumin peptide OVA (257-264) class I MHC molecule
AV-9320-1	Ovalbumin peptide OVA (323-339) MHC class II peptide
AV-9325-10	Dinitrophenyl (DNP)-KLH protein Conjugate
AV-9335-10	Dinitrophenyl (DNP)-Ovalbumin (OVA) protein Conjugate
AV-9340-1	Dinitrophenyl (DNP)-Lipopolysaccharide (LPS) Conjugate

Complete list is available at:
http://4adi.com/objects/catalog/product/extras/Vaccine_Adjuvants_flr.pdf

AV-9320-1 140518P

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com