




Life Technologies™

Your Molecular & Cell Technology Partner



Buffer Pack, DAS or TAS ELISA, Alkphos, No ECI Buffer



Intended Use:

This buffer pack does not include ECI buffer but contains all other buffers and substrate used in most DAS and TAS alkphos labeled ELISA systems.

Includes:

- Carbonate coating buffer, 10X concentrate
- PNP substrate buffer, 5X concentrate
- PNP substrate tablets, 5 mg each
- PBST buffer, powder
- General extract buffer
- Tween 20
- [User Guide](#)

Option
Buffer Pack, DAS or TAS ELISA without ECI buffer, Alkphos, 500 tests ACC 00333/0500
Buffer Pack, DAS or TAS ELISA without ECI buffer, Alkphos, 1000 tests ACC 00333/1000
Buffer Pack, DAS or TAS ELISA without ECI buffer, Alkphos, 5000 tests ACC 00333/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com