

# Beaker, Graduated, 4000 mL

Product ID B199

## Description

Graduated, Polypropylene, 1000 mL graduations.

## Product Information

Autoclavability	Up to 121 °C
Capacity	4000 mL
Graduations	1000 mL
Material	Polypropylene
Product Number	B199

## Shipping and Storage

Storage Temp.	Room Temperature
Tariff Code	3923.30.0000

**PhytoTechnology Laboratories®**

14610 W 106th St. Lenexa, KS 66215

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)