

# Beaker, 800 mL, Tri-Corner

Product ID B223

## Description

Three corner spouts for easy pouring.

## Product Information

Autoclavability	Up to 121 °C
Capacity	800 mL
Graduations	50 mL
Material	Polypropylene
Product Number	B223

## Shipping and Storage

Storage Temp.	Room Temperature
Tariff Code	3926.90.9910

**PhytoTechnology Laboratories®**

14610 W 106th St. Lenexa, KS 66215

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)