

Beaker, 1000 mL, Tri-Corner

Product ID B225

Description

Three corner spouts for easy pouring.

Product Information

Autoclavability	Up to 121 °C
Capacity	1000 mL
Graduations	50 mL
Material	Polypropylene
Product Number	B225

Shipping and Storage

Storage Temp.	Room Temperature
Tariff Code	3926.90.9910

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com