

Beaker, Graduated, 1000 mL

Product ID B931

Description

Polymethyl pentene - Not autoclavable.

Product Information

Autoclavability	Up to 121 °C
Capacity	1000 mL
Graduations	50 mL
Material	Polymethyl Pentene
Product Number	B931

Shipping and Storage

Storage Temp.	Room Temperature
Tariff Code	3923.30.0000

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com