

Beaker, Graduated, 5000 mL

Product ID B960

Description

Polymethyl pentene - Not autoclavable.

Product Information

| | |
|-----------------|--------------------|
| Autoclavability | Up to 121 °C |
| Capacity | 5000 mL |
| Material | Polymethyl Pentene |
| Readability | 500 mL |
| Product Number | B960 |

Shipping and Storage

| | |
|---------------|------------------|
| Storage Temp. | Room Temperature |
| Tariff Code | 3923.30.0000 |

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com