

Wash Bottle, Isopropyl Alcohol

Product ID B974

Description

Marked bottle for isopropyl alcohol only
Polyethylene
500 mL capacity.

Product Information

Autoclavability	Not autoclavable
Capacity	500 mL
Use With	Isopropyl Alcohol
Material	Polyethylene
Product Number	B974

Shipping and Storage

Storage Temp.	Room Temperature
Tariff Code	3923.90.0080

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com