

Wash Bottle, Distilled Water

Product ID B987

Description

Marked bottle for distilled water only
Polyethylene
500 mL capacity.

Product Information

Autoclavability	Not autoclavable
Capacity	500 mL
Use With	Distilled Water
Material	Polyethylene
Product Number	B987

Shipping and Storage

Storage Temp.	Room Temperature
Tariff Code	3923.90.0080

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com