

Bacterial Reagent Set for *Xanthomonas axonopodis* pv. *begoniae* (Xab)

This ELISA test is a qualitative serological assay for the detection of *Xanthomonas axonopodis* pv. *begoniae* (Xab) in bacterial cultures and plant vascular tissue streams. Xab is a member of the *Xanthomonas* genus known for having 27 plant disease causing species able to infect 400 different plant species.

If you would prefer to add a field-deployable assay to your testing toolkit, please check out Agdia's [Xan Immunoseras](#)® (will detect genus and not differentiate species).

Included:

- Detection antibody
- Alkaline phosphatase enzyme conjugate
- 96-well microtiter plates, solid* (*strip in single plate set)
- [User Guide](#)

Option

Bacterial Reagent Set, *Xanthomonas axonopodis* pv. *begoniae* (Xab), 96 testwells

BRA 23700/0096

Bacterial Reagent Set, *Xanthomonas axonopodis* pv. *begoniae* (Xab), 500 testwells

BRA 23700/0500

Bacterial Reagent Set, *Xanthomonas axonopodis* pv. *begoniae* (Xab), 1000 testwells

BRA 23700/1000

Bacterial Reagent Set, *Xanthomonas axonopodis* pv. *begoniae* (Xab), 5000 testwells

BRA 23700/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com