

## Closure, Round, Screw Cap, 70 mm, Vented

Product ID: **C170**

### Introduction

---

**Description:** Natural polypropylene closure for 70 mm diameter screw-cap culture vessels (glass and plastic). Vented (0.5 micron pore size) and autoclavable

#### Product Information

Autoclavability	Up to 121 °C
Diameter	70 mm
Use With	C583, C590
Material	Polypropylene
Storage Temp	Room Temperature
UPC / SKU	C170



#### India Contact:



#### Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, opposite M2K Pitampura, Delhi-110034 ( India)  
Tel # +91-11-4220-8000; 4220-8111; 4220-8222 Fax# +91-11-4220-8444, Mobile# +91-98105-21400,  
Toll Free # 1800-120-2434

Email - [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) | [customerservice@atzlabs.com](mailto:customerservice@atzlabs.com)

ISO Certifications: [ISO 9001:2008 \(QMS\)](#) | [ISO 10002:2014](#)

Life Technologies, Lifetechindia.com, Lifetechindia.co.in , Lifetechnologiesindia.com , Lifetechnologiesindia.co.in, LifeTech, LifeTech (triangular arc), Auprep, Auprep.com, AuPreP Oligos Gelrose, Novascript, neofect**AMINE 3000**, genfect**AMINE 4000**, ATZ labs, ATZ Touch, ATZ Dew, ATZ Fabrich, ATZ MobiFriend, ATZLabs.com are the key trademarks/domains of Life Technologies India Private Limited and its associated entities.

Important: We have not granted licence or permission to any person or entity to use our above trade marks / domains  
Any person / entity found to be using our marks shall be liable for civil / criminal action.

**All Products shown on our website are for Laboratory research applications unless otherwise stated**

© 2014. Life Technologies (India) Pvt. Ltd. All rights reserved.  
design by [www.nextconcept.in](http://www.nextconcept.in)