

Culture Vessel, Glass, 70 mm, 8 oz (240 mL)

Product ID c1892

Description

Culture Vessel, Glass, 70 mm, 8 oz (240 mL) C1892 glass culture vessels are compatible with our 70 mm closures (included). This product replaces C583 (see dimensions and specifications in the Product Information tab)..

Product Information

Autoclavability	Up to 121 °C
Diameter	66 mm (base), 72 mm (top)
Use With	C566 (included), C170 (sold separately)
Material	Glass
Opening Diameter	70 mm
Product Number	C583
Height	90 mm

Shipping and Storage

Tariff Code	7010.90.5057
-------------	--------------

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com