

## Culture Vessel, SteriCon™-13

Product ID C1932


### Description

SteriCons are square, sterilized culture vessels. They are ready to be used straight out of the box. Lids and container bases are packaged separately in their own sleeves.

**Note:**

*NOT autoclavable*

### Product Information



Autoclavability	Not autoclavable
Volume	946 mL
Material	Polystyrene
Product Number	C1932
Length	101.6 mm (top)
Width	101.6 mm (top)
Height	130 mm

**India Contact:**



**Life Technologies (India) Pvt Ltd.**

306, Agarwal City Mall, opposite M2K Pitampura, Delhi-110034 ( India)  
Tel # +91-11-4220-8000; 4220-8111; 4220-8222 Fax# +91-11-4220-8444, Mobile# +91-98105-21400,  
Toll Free # 1800-120-2434  
Email - [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) | [customerservice@atzlabs.com](mailto:customerservice@atzlabs.com)  
ISO Certifications: [ISO 9001:2008 \(QMS\)](#) | [ISO 10002:2014](#)