

Culture Tube, PhytoTech® Brand, 25 mm x 95 mm Round-Bottom, Polypropylene

Product ID C2039

Description

Polypropylene culture tubes are autoclavable (reusable) and very durable. No need to worry about broken glass when manipulating closures or dropping them on your work surface.

Our polypropylene tubes are compatible with all culture tube closures sold by *PhytoTechnology Laboratories*. (***sold separately***)

We recommend shorter culture tube racks, like our stainless steel rack or 12-place no-wire rack, for use with shorter culture tubes like this one.

Product Information

Autoclavability	Up to 121 °C
Use With	C1800, C1805, C069, C163, C945 (all sold separately)
Volume	30 mL
Material	Polypropylene
Opening Diameter	25 mm
Product Number	C2040
Height	150 mm