

Culture Tube, PhytoTech® Brand, 25 mm x 150 mm Round- Bottom, Polypropylene

Product ID c2040

Description

Polypropylene culture tubes are autoclavable (reusable) and very durable. No need to worry about broken glass when manipulating closures or dropping them on your work surface.

Our polypropylene tubes are compatible with all culture tube closures and culture tube racks sold by *PhytoTechnology Laboratories*. **(all sold separately)**

)
..

Product Information

Autoclavability	Up to 121 °C
Use With	C1800, C1805, C069, C163, C945 (all sold separately)
Volume	55 mL
Material	Polypropylene
Opening Diameter	25 mm
Product Number	C2040
Height	150 mm

Shipping and Storage

Tariff Code	3926.90.9900
-------------	--------------

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com