

Culture Vessel, PTL-100™

Product ID c2100

Description

The PTL-100 culture vessel is a square box-style culture vessel that is compact enough for optimal storage, yet large enough to allow for multiple cultures to grow in your gel. PTL-100 is also graduated with 25 mL and 50 mL graduations on one side so you have a guide for how much media you are pouring into your vessel. PTL-100's are molded from autoclavable high clarity polycarbonate plastic which allows for lightweight transportation of your cultures, excellent light transmission, and high durability (virtually unbreakable). PTL-100's are packaged with PTL-100C closures.

Please note: *The container and lid is very tight fitting at first, but loosens up after the first autoclave cycle*

..

Product Information

| | |
|-----------------|------------------|
| Autoclavability | Up to 121 °C |
| Use With | C2101 (included) |
| Volume | 372 mL |
| Material | Polycarbonate |
| Length | 75 mm |
| Width | 75 mm |
| Height | 98 mm |

Shipping and Storage

| | |
|-------------|--------------|
| Tariff Code | 3926.90.9910 |
|-------------|--------------|

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com