

Closure, PTL-100C™

Product ID c2101

Description

Closure for PTL-100 culture vessels

Product Information

Autoclavability	Up to 121 °C
Use With	C2100 (sold separately)
Material	Polypropylene
Length	75 mm
Width	75 mm

Shipping and Storage

Tariff Code	3926.90.9910
-------------	--------------

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com