

# Closure, PTL-100C™, Vented

Product ID c2110

## Description

Closure for PTL-100 culture vessels. Vented (19 mm vent; 0.5 micron pore size).

## Product Information

|                 |                         |
|-----------------|-------------------------|
| Autoclavability | Up to 121 °C            |
| Use With        | C2100 (sold separately) |
| Material        | Polypropylene           |
| Length          | 75 mm                   |
| Width           | 75 mm                   |

## Shipping and Storage

|             |              |
|-------------|--------------|
| Tariff Code | 3926.90.9910 |
|-------------|--------------|

**PhytoTechnology Laboratories®**

14610 W 106th St. Lenexa, KS 66215

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)