

Culture Vessel, SteriCon™-8

Product ID: **C2118**

Introduction

Description: SteriCons are square, sterilized culture vessels. They are ready to be used straight out of the box. Lids and container bases are packaged separately in their own sleeves.

Note: NOT autoclavable

Product Information

Autoclavability	Not autoclavable
Volume	473 mL
Material	Polystyrene
Storage Temp	Room Temperature
UPC / SKU	C2118
Length	101.6 mm (top)
Width	101.6 mm (top)
Height	74.9 mm



India Contact:



Life Technologies (India) Pvt Ltd.

306, Agarwal City Mall, opposite M2K Pitampura, Delhi-110034 (India)

Tel # +91-11-4220-8000; 4220-8111; 4220-8222 Fax# +91-11-4220-8444, Mobile# +91-98105-21400,

Toll Free # 1800-120-2434

Email - customerservice@lifetechindia.com | customerservice@atzlabs.com

ISO Certifications: [ISO 9001:2008 \(QMS\)](#) | [ISO 10002:2014](#)

