

Product Specification Sheet
Human Cancer Antigen 19-9 (CA19-9) Antibody

Cat. #	CA1991-M
Lot #	148
Description	Mouse Monoclonal Anti-Human CA19-9 IgG # 1
Antigen	Highly purified CA19-9 protein, Human
Host Animal	Mouse
Isotype	IgG1
Form	Protein A/G affinity purified from ascites fluid supplied in PBS, pH 7.4 and no preservative
Antibody concentration	2.6 mg/ml
Storage	Short term at 2-8oC; Long term at -20oC
Stability	~1 year or more at -20oC; Avoid repeated freeze and thaw.
Specificity and crossreactivity	Specific for sialylated Lewis ^a (SLe ^a)
Suggested Usage	Suitable for ELISA as coating antibody and tracer antibody; concentrations and conditions must be optimized for the assay
MSDS info	Natural protein; MSDS not applicable

CA1991-M

71220A