

## Coating Antibody for Strawberry mild yellow edge virus (SMYEV)



**Test Format:** Triple Antibody Sandwich (TAS) ELISA

### Intended Use:

This capture antibody is intended for use with its corresponding detection antibodies (ECA).

Test performance using other detection antibodies has not been evaluated and is not guaranteed.

### Includes:

- Capture Antibody
- 96-well microtiter plates, solid

### Option

Coating antibody for Strawberry mild yellow edge virus (SMYEV), 500 testwells

CAB 14500/0500

Coating antibody for Strawberry mild yellow edge virus (SMYEV), 1000 testwells

CAB 14500/1000

Coating antibody for Strawberry mild yellow edge virus (SMYEV), 5000 testwells

CAB 14500/5000

### India Contact:

#### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)