

Coating Antibody for *Xanthomonas hortorum* pv. *pelargonii* (Xhp)



Test Format: Double Antibody Sandwich (DAS) ELISA
Synonym: *Xanthomonas campestris* pv. *pelargonii* (Xcp)

Intended Use:

This capture antibody is intended for use with its corresponding enzyme conjugate (ECA).

Test performance using other enzyme conjugate(s) has not been evaluated and is not guaranteed.

Includes:

- Capture Antibody
- 96-well microtiter plates, solid

Option
Coating antibody for <i>Xanthomonas hortorum</i> pv. <i>pelargonii</i> (Xhp), 500 testwells CAB 32503/0500
Coating antibody for <i>Xanthomonas hortorum</i> pv. <i>pelargonii</i> (Xhp), 1000 testwells CAB 32503/1000
Coating antibody for <i>Xanthomonas hortorum</i> pv. <i>pelargonii</i> (Xhp), 5000 testwells CAB 32503/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com