

Coating Antibody for *Xanthomonas axonopodis* pv. *phaseoli* (Xaph)



Test Format: Double Antibody Sandwich (DAS) ELISA

Intended Use:

This capture antibody is intended for use with its corresponding enzyme conjugate (ECA).

Test performance using other enzyme conjugate(s) has not been evaluated and is not guaranteed.

Includes:

- Capture Antibody
- 96-well microtiter plates, solid

Option

Coating antibody for *Xanthomonas axonopodis* pv. *phaseoli* (Xaph), 500 testwells
CAB 36101/0500

Coating antibody for *Xanthomonas axonopodis* pv. *phaseoli* (Xaph), 1000 testwells
CAB 36101/1000

Coating antibody for *Xanthomonas axonopodis* pv. *phaseoli* (Xaph), 5000 testwells
CAB 36101/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com