

Coating Antibody for *Pseudomonas savastanoi* pv. *phaseolicola* (Psph)



Test Format: Double Antibody Sandwich (DAS) ELISA
Synonym: *Pseudomonas syringae* pv. *phaseolicola*

Intended Use:

This capture antibody is intended for use with its corresponding enzyme conjugate (ECA).

Test performance using other enzyme conjugate(s) has not been evaluated and is not guaranteed.

Includes:

- Capture Antibody
- 96-well microtiter plates, solid

Option

Coating antibody for *Pseudomonas savastanoi* pv. *phaseolicola* (Psph), 500 testwells
CAB 53300/0500

Coating antibody for *Pseudomonas savastanoi* pv. *phaseolicola* (Psph), 1000 testwells
CAB 53300/1000

Coating antibody for *Pseudomonas savastanoi* pv. *phaseolicola* (Psph), 5000 testwells
CAB 53300/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com