

Detection Antibody for *Xanthomonas oryzae* pv. *oryzae* (Xoo)



Test Format: Antigen coated plate (ACP) ELISA
Synonym: *Xanthomonas campestris* pv. *oryzae*

Intended Use:

This detection antibody is intended for use with its corresponding enzyme conjugate (ECA).

Test performance using other enzyme conjugate(s) has not been evaluated and is not guaranteed.

Includes:

- Detection Antibody
- 96-well microtiter plates, solid

Option

Detection antibody for *Xanthomonas oryzae* pv. *oryzae* (Xoo), 500 testwells
CAB 85000/0500

Detection antibody for *Xanthomonas oryzae* pv. *oryzae* (Xoo), 1000 testwells
CAB 85000/1000

Detection antibody for *Xanthomonas oryzae* pv. *oryzae* (Xoo), 5000 testwells
CAB 85000/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com