

Catalog#	CD03-D5-100	Size	100 tests
Catalog#	CD03-D5-25	Size	25 tests

Mouse Anti Human CD3 FITC/CD8 PE

Clone Number: HIT3b/HIT8a

Isotype: IgG1/IgG1

Volume/Quantity: 50T

Product Form: Dual color combination consisting of FITC conjugated and R-PE conjugated monoclonal antibodies mixed in optimal ratio.

USAGE

The reagent is tested for flow cytometric analysis using 20µl/10⁶ cells or 100µl peripheral blood cells.

STORAGE

Storage at 4°C, should not be frozen and avoid prolonged exposure to light.

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com