

Catalog# CD03-T2-100 Size 50 tests
Catalog# CD03-T2-25 Size 25 tests

Mouse Anti Human CD3 FITC/CD4 PE/CD45 PE-Cy5

Clone Number: HIT3b/HIT4a/HI30

Isotype: IgG1/IgG2b/IgG1

Product Form: Triple color combination consisting of FITC, R-PE and PE-Cy5 conjugated monoclonal antibodies mixed in optimal ratio.

USAGE

The reagent is tested for flow cytometric analysis using 20µl/10⁶ cells or 100µl peripheral blood cells.

STORAGE

Storage at 4°C, should not be frozen and avoid prolonged exposure to light.

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com