

<b>Catalog#</b>	<b>CD04-T1-100</b>	<b>Size</b>	<b>50 tests</b>
<b>Catalog#</b>	<b>CD04-T1-25</b>	<b>Size</b>	<b>25 tests</b>

---

**Mouse Anti Human CD4 FITC/CD8 PE/CD3 PE-Cy5**

**Clone Number:** HIT4a/HIT8a/HIT3b

**Isotype:** IgG2b/IgG1/IgG1

**Volume/Quantity:** 25T

**Product Form:** Triple color combination consisting of FITC, R-PE and PE-Cy5 conjugated monoclonal antibodies mixed in optimal ratio.

**USAGE**

The reagent is tested for flow cytometric analysis using 20µl/10<sup>6</sup> cells or 100µl peripheral blood cells.

**STORAGE**

Storage at 4oC, should not be frozen and avoid prolonged exposure to light.

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)