
Mouse Anti Human CD5 FITC/CD19 PE

Catalog # CD05-D1-100 **Size** 100 tests
Catalog # CD05-D1-25 **Size** 25 tests

Clone Number: HISM2/HIB19

Isotype: IgG1/IgG1

Volume/Quantity: 25T

Product Form: Dual color combination consisting of FITC conjugated and R-PE conjugated monoclonal antibodies mixed in optimal ratio.

USAGE

The reagent is tested for flow cytometric analysis using 20µl/10⁶ cells or 100µl peripheral blood cells.

STORAGE

Storage at 4oC, should not be frozen and avoid prolonged exposure to light.

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com