

---

## Mouse Anti-Human CD5-PE-Cy5-conjugate

**Catalog #** CD05PC-100 Size 100 tests  
**Catalog #** CD05PC-25 Size 25 tests

**Clone Number:** HISM2/HIB19

**Isotype:** IgG1/IgG1

**Volume/Quantity:** 25T

**Product Form:** Dual color combination consisting of FITC conjugated and R-PE conjugated monoclonal antibodies mixed in optimal ratio.

### USAGE

The reagent is tested for flow cytometric analysis using 20µl/10<sup>6</sup> cells or 100µl peripheral blood cells.

### STORAGE

Storage at 4oC, should not be frozen and avoid prolonged exposure to light.

India Contact:

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)