
Mouse Anti Human CD8 FITC/CD2 PE

Catalog # **CD08-D1-100** **Size# 100 tests**
Catalog # **CD08-D1-25** **Size# 25 tests**

Clone Number: HIT8a/HIT11

Isotype: IgG1/IgG1

Product Form: Dual color combination consisting of FITC conjugated and R-PE conjugated monoclonal antibodies mixed in optimal ratio.

USAGE

The reagent is tested for flow cytometric analysis using 20µl/10⁶ cells or 100µl peripheral blood cells.

STORAGE

Storage at 4°C, should not be frozen and avoid prolonged exposure to light.

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com