

---

## Mouse Anti Human CD19 FITC/CD2 PE

<b>Catalog #</b>	<b>CD19-D1-100</b>	<b>Size</b>	<b>100 tests</b>
<b>Catalog #</b>	<b>CD19-D1-25</b>	<b>Size</b>	<b>25 tests</b>

**Clone Number:** HIB19/HIT11

**Isotype:** IgG1/IgG1

**Product Form:** Dual color combination consisting of FITC conjugated and R-PE conjugated monoclonal antibodies mixed in optimal ratio.

### USAGE

The reagent is tested for flow cytometric analysis using 20µl/10<sup>6</sup> cells or 100µl peripheral blood cells.

### STORAGE

Storage at 4°C, should not be frozen and avoid prolonged exposure to light.