
Mouse Anti Human CD19 FITC/CD3 PE

Category	CD19-D2-100	Size	100 tests
Category	CD19-D2-25	Size	25 tests

Clone Number: HIB19/HIT3b
Isotype: IgG1/IgG1

Product Form: Dual color combination consisting of FITC conjugated and R-PE conjugated monoclonal antibodies mixed in optimal ratio.

USAGE

The reagent is tested for flow cytometric analysis using 20µl/10⁶ cells or 100µl peripheral blood cells.

STORAGE

Storage at 4°C, should not be frozen and avoid prolonged exposure to light.